## Polyvinylpyrrolidone (PVP) 7% (w/v) in HEPES-HTF

**For laboratory procedures only; other uses must be qualified by the end user.**

Product Description	REF Number	Unit Size
PVP 7% (w/v) in HEPES-HTF	ART-4005-A	6 x 0.5 mL

#### INTENDED USE

Polyvinylpyrrolidone (PVP) of average molecular weight 360,000 has been used for the ICSI procedure to decrease the motility of spermatozoa, prevent their sticking to the injection pipette during the procedure, and to give more control over the flow of fluid in the ICSI needle (Palermo et al. Fertil Steril. 1993;59:826-835). These features are due primarily to the increased viscosity of the PVP solution.

#### PRODUCT DESCRIPTION

PVP at 7% concentrations yield a viscosity that is suitable for sperm function as well as improved visibility to detect sperm.

This product contains 10 mg/L of gentamicin, an aminoglycoside antibiotic.

#### COMPONENTS

PVP with an average molecular weight of 360,000 in HEPES-HTF contains 5 mg/mL Human Serum Albumin, and is filtered through a sterile 0.2-µ filter.

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#### PRECAUTIONS AND WARNINGS

Do not use medium that shows evidence of particulate matter, cloudiness, or is not rose colored.

**Caution:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when testing for antibodies to HIV-1/HIV-2, HCV and non-reactive for HBsAg, HCV RNA and HIV-1 RNA. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents. Donors of the source material have been screened for Creutzfeldt-Jacobs disease (CJD). Based on effective donor screening and product manufacturing processes, it carries an extremely remote risk for transmission of viral diseases. A theoretical risk for transmission of CJD is also considered extremely remote. No cases of transmission of viral diseases or CJD have ever been identified for albumin.

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

The PVP 7% (ART-4005) component in this kit contains 5 mg/mL human serum albumin.

**Single use:** To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

**Note:** Embryo is considered a general term. More precisely, SAGE considers the period of time initiating when a single diploid cell results from the fusion of male and female genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation (typically day 5-6) ending with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

**Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

This product contains the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

#### QUALITY ASSURANCE

One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

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A certificate of analysis is available for this product.

#### DIRECTIONS FOR USE

A washed sperm preparation can be diluted to 1 to 2 x 10<sup>6</sup> motile sperm in 0.5 to 1.0 mL of Quinn's® Sperm Washing Medium (REF #1005/1006). This sample is spun at high speed (approx. 1800 x g for 5 minutes immediately prior to use). The sperm pellet is suspended in the same supernatant in the tube. The sperm preparation for ICSI is then made by adding 1 to 2 µL of the sperm suspension to a 5 to 10 µL drop of PVP solution under oil (REF #4008) in the chamber used for ICSI. Volumes of sperm preparation and PVP solution, as well as sperm concentrations, may have to be modified to give optimal final sperm density for ICSI.

**Each laboratory should make its own determination of which medium to use for each particular procedure.**

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog.

#### STORAGE INSTRUCTIONS AND STABILITY

Store unopened containers refrigerated at 2 °C to 8 °C. Warm to incubator (37 °C) temperature prior to use. Exposure of the medium to CO<sub>2</sub> should be minimized to avoid pH levels of 7.0 or less. Do not freeze or expose to temperatures greater than 39 °C. The product is stable until the expiration date shown on the label.

- Remove desired volume of product using aseptic procedures.
- Once removed, do not return any volume of product to the original container.
- Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

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ART-1005/1006 Quinn's® Sperm Washing Medium

ART-4008P Oil for Tissue Culture

SAGE In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE Support Line.

Quinn's® is a registered trademark of CooperSurgical, Inc.

**Call the SAGE SUPPORT LINE: In the U.S. : (800) 243-2974 International: (203) 601-9818**

**Each laboratory should make its own determination of which medium to use for each particular procedure.**

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog.

#### STORAGE INSTRUCTIONS AND STABILITY

Store unopened containers refrigerated at 2 °C to 8 °C. Warm to incubator (37 °C) temperature prior to use. Exposure of the medium to CO<sub>2</sub> should be minimized to avoid pH levels of 7.0 or less. Do not freeze or expose to temperatures greater than 39 °C. The product is stable until the expiration date shown on the label.

- Remove desired volume of product using aseptic procedures.
- Once removed, do not return any volume of product to the original container.
- Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

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#### EXPLANATION OF SYMBOLS

	Catalog Number
	Batch Number
	Use By (year, month, day)
	Do Not Reuse
	Temperature Limitation
	Aseptic Technique Sterilization Membrane Filtered (SAL 10 <sup>-3</sup> )
	<b>ATTENTION:</b> See instructions for use.
	Authorized Representative in the European Community.
	Product conforms to the Medical Device Directive 93/42/EEC
	Manufacturer
	U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

	SAGE In Vitro Fertilization <i>a CooperSurgical Company</i>
	SAGE In Vitro Fertilization, Inc. a CooperSurgical Company 95 Corporate Drive Trumbull, CT 06611 USA
	Leisegang Feinmechanik GmbH Leibnizstraße 32 D-10625, Berlin GERMANY www.coopersurgical.com
	SAGE In Vitro Fertilization 1979 East Locust Street Pasadena, CA 91107 USA
	SPN-00595 • Rev. D 2/11

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	SAGE In Vitro Fertilization <i>a CooperSurgical Company</i>
	SAGE In Vitro Fertilization, Inc. a CooperSurgical Company 95 Corporate Drive Trumbull, CT 06611 USA
	Leisegang Feinmechanik GmbH Leibnizstraße 32 D-10625, Berlin GERMANY www.coopersurgical.com
	SAGE In Vitro Fertilization 1979 East Locust Street Pasadena, CA 91107 USA
	SPN-00595 • Rev. D 2/11

## Italiano

## Polyvinylpyrrolidone (PVP) 7% (w/v) in HEPES-HTF (Polivinilpirrolidone (PVP) 7% (w/v) in HEPES-HTF)

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