# FertiCult™ IVF medium

Cell culture medium for in vitro culture of human embryos, during the first 48 hours in culture

Document reference: FP09 I07 R01 B.4, Update: 24/OCT/2013

## GENERAL INFORMATION

FertiCult IVF medium is a ready to use formulation for the *in vitro* culture of mammalian embryos. It is designed for short term culture only (up to 48 hours in culture). The medium is complete and needs no further additives. If preferred 10 % patient serum v/v can be added.

FertiCult IVF medium is suitable for micro-droplets under oil or equally suited to culture in Falcon 3037 organ culture dishes (also 4-well Nunc culture dishes).

As with all IVF media, FertiCult IVF medium has to be **pre-incubated in CO<sub>2</sub> incubator with 5-6% CO<sub>2</sub> for 24 hours before use** (with lid opened).

## MATERIAL INCLUDED WITH THE KIT

Product code	Product description
FECU020	5 x 20ml FertiCult IVF medium
FECU050	5 x 50mL FertiCult IVF medium
FECU100	3 x 100mL FertiCult IVF medium
FECU020PHR	5 x 20ml FertiCult IVF medium w/ phenol red
FECU050PHR	5 x 50mL FertiCult IVF medium w/ phenol red
FECU100PHR	3 x 100mL FertiCult IVF medium w/ phenol red
FECU020PHR_G	5 x 20ml FertiCult IVF medium w/ phenol red and Gentamicin
FECU050PHR_G	5 x 50mL FertiCult IVF medium w/ phenol red and Gentamicin
FECU100PHR_G	3 x 100mL FertiCult IVF medium w/ phenol red and Gentamicin

## MATERIAL NOT INCLUDED WITH THE KIT

- Incubator at 37°C (5-6%CO<sub>2</sub>)
- Petri dishes (e.g. Falcon 3037)
- Mineral oil (e.g. FertiCult Mineral Oil)
- Laminar flow bench (ISO5 environment)
- Microscope
- Test tubes

## PRODUCT SPECIFICATIONS

- Chemical composition
- pH: between 7.20 - 7.50 (37°C – 5-6% CO<sub>2</sub>)
- Osmolality: 270 - 290 mOsm/kg
- Sterility: sterile (SAL 10<sup>-3</sup>)
- Endotoxins: < 0.25 EU/ml
- Mouse Embryo Assay (blastocysts after 96h culture): ≥ 80% (48h exposure from zygote stage)
- Use of Ph Eur or USP grade raw materials if applicable
- Certificate of analysis and MSDS are available upon request

## PRE-USE CHECKS

- Do not use the product if it becomes discoloured (if medium contains phenol red), cloudy, or shows any evidence of microbial contamination
- Do not use the product if seal of the container is opened or defect when the product is delivered

## STORAGE AND CONSERVATION

- Store between 2-25°C, once opened store between 2-8°C
- Do not freeze before use
- Keep away from (sun)light
- After opening the container do not use the product longer than 7 days
- Do not use after expiry date

## WARNINGS AND PRECAUTION

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be

totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Always wear protective clothing when handling specimens.

Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination, even when FertiCult IVF medium contains Gentamicin.

FertiCult IVF medium with phenol red and gentamicin contains the antibiotic Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

## METHOD

### PREPARATIONS FOR USE OF MICRO-DROPLETS

For micro-droplets, between 100-250µl of FertiCult IVF medium may be dispensed around the culture dish, up to 8 per 60mm dish.

The dish is then filled with 5ml of pre-washed and pre-equilibrated light mineral oil (nontoxic and preferably embryo tested, e.g. FertiCult Mineral Oil).

One oocyte is usually placed with about 10 000 sperm cells per micro-droplet, so if each patient during the next day is expecting 16 oocytes, at least two dishes per patient should be prepared. The dishes are then placed in the incubator (usually non-humidified) to equilibrate overnight at 37°C and under an atmosphere of 5-6% CO<sub>2</sub> in air.

**Overnight equilibration is highly recommended.**

### PREPARATIONS FOR USE IN OPEN SYSTEMS

In open systems such as with the Falcon or Nunc dishes, about 1ml of medium is placed in each well.

A further 3ml is placed in the reservoir surrounding the wells. This helps to maintain humidity as well as providing medium with which to wash the oocytes and embryos.

Up to five oocytes per dish (or well if using Nunc dishes) may be cultured together, usually with about 100 000 sperm cells per well.

So if each patient during the next day is expecting 15 oocytes at oocyte collection, at least six dishes per patient should be prepared (three for initial holding and cleaning of oocytes during collection and then three dishes for overnight culture).

When using Nunc dishes, two would be enough (since there are 4 wells per dish, the resulting 8 wells would be more than enough).

The dishes are then placed in the incubator to equilibrate overnight at 37°C and under an atmosphere of 5-6% CO<sub>2</sub> in air.

As with the micro-droplets method described above, **equilibration overnight is highly recommended.**

In the open system, paraffin oil is not necessary, but sometimes a 1ml layer of mineral oil may be added over the medium in the inner well. When using such an open system the incubator must be humidified.

## CULTURE

After equilibration overnight, the prepared dishes are ready for use.

During oocyte collection, oocytes are identified in the follicular aspirates and then washed before placing into the wells or droplets in the dish (as appropriate). Once oocyte collection is over, oocytes will be checked and may be rewashed, before then being placed into fresh droplets/wells for overnight culture.

Generally, within 6 hours of oocyte collection, oocytes will be inseminated with prepared sperm and then left overnight in the incubator.

The procedures described above are for the first 24 hours of culture.

Once embryos have been identified by the presence of two pro-nuclei, they are usually placed into fresh dishes containing fresh medium.

In the past, where people have used cord serum or patient serum, the initial medium contains 10% and the changeover (after fertilisation has been confirmed) 15%.

When using FertiCult IVF medium no extra protein is required at changeover. Dishes should be prepared as described above, equilibrated overnight, after which they will be ready to receive the newly fertilised embryos.

**NB.** Test tube culture is now rare, but for those still using this method, place 1ml of FertiCult IVF medium in each tube and then follow the same procedure as for Falcon or Nunc dishes.

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