## SpermFreeze™ SSP

Doc. ref. FP09 I11 R01-SSP A.4, Update: 15/02/2013  
SpermFreeze SSP is sterilized by sterile filtration

### INTRODUCTION

SpermFreeze SSP is a HEPES buffered freezing medium for use with human sperm<sup>1,2</sup>.

SpermFreeze SSP contains 0.4 % human serum albumin and 26.76% glycerol to protect the sperm from damage due to the freezing procedure<sup>3</sup>.

### SPERMFREEZE SSP AND SPERM PREPARATION

SpermFreeze SSP can be used in combination with Sil-Select Plus™<sup>4,5</sup>.

#### Before freezing

In case of very low sperm concentrations it is advisable to concentrate the sperm before freezing. This may increase sperm quality after thawing and will reduce the number of straws to be frozen.

#### After thawing

Use sperm preparation techniques after thawing to eliminate dead sperm cells and debris. Dilute the concentrated sperm in washing medium or any other medium you would normally use.

### SPECIFICATIONS AND QUALITY CONTROL

- pH: 7.20 - 7.60
- Sterility: sterile
- Endotoxin: < 0.25EU/ml
- Sperm survival test: ≥ 80% survival after 45 minute exposure of untreated semen to the test medium
- Not MEA tested
- A certificate of analysis is available upon request.

### METHOD

#### Freezing

1. Allow the semen to liquefy at room temperature for 30 minutes.
2. Mix 3 parts sperm with 1 part SpermFreeze SSP. Add the medium in drops while gently swirling.  
**Caution:** To avoid cold-shock, make sure SpermFreeze SSP is at room temperature.
3. Leave the mixture for 10 minutes at room temperature for equilibration.
4. Suck the sample/medium mixture into the freezing straws, leaving approximately 1.5cm of air at the end of the straw.
5. Seal the straws.
6. Dry off individually with a linen free cloth.
7. Shake to move the air-bubble to the centre of the straw.
8. Freeze vertically for 15 minutes, just above the level of the liquid nitrogen.
9. Store in liquid nitrogen.

#### Thawing

1. Remove as many straws as required from the liquid nitrogen.
2. Place the straws in tap water for 5 minutes.
3. Cut off the end of the straw, place the open end inside a container (e.g. a test tube) and tap the straw against the side of the container to allow complete evacuation of the mixture.
4. Dilute the concentrated sperm in a suitable insemination medium (at least 3 ml per 0.5 ml semen) and mix thoroughly
5. Centrifuge during 15 minutes at 300-350g
6. Resuspend pellet in a suitable insemination medium (e.g. FertiCult Flushing medium)

### MATERIAL INCLUDED IN THE KIT

#### SSP001

- 5x 1ml of SpermFreeze with 26.76% glycerol

### MATERIAL NOT INCLUDED IN THE KIT

- Plastic freezing straws

### STORAGE AND CONSERVATION

Store SpermFreeze SSP between 2-8°C upon receipt.

Keep from light.

### WARNINGS AND PRECAUTIONS

All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.

**Caution:** All blood products should be treated as potentially infectious. Source materials from which this product was derived was found negative when tested for antibodies to HIV, HBc, HCV, and HTLV I/II and non-reactive for HBsAg, HCV RNA, and HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents.

### REFERENCES

- 1 Mahadevan M, Trounson AD. Effect of cryoprotective media and dilution methods on the preservation of human spermatozoa. *Andrologia*, 1983; 15: 355-66.
- 2 Mahadevan M, Trounson AD, Leeton JF. Successful use of human semen cryobanking for in vitro fertilization, *Fertil Steril*, 1983; 15: 355-66.
- 3 Brotherton J. Cryopreservation of human semen. *Archives of Andrology*, 1990; 25: 181-95.
- 4 Kobayashi T, Kaneko S, Hara I, Park YJ, et al. Concentrating human sperm before cryopreservation. *Andrologia*, 1991; 23: 25-8.
- 5 Graczykowski JW, Siegel MS. Influence of sperm processing on the fertilizing capacity and recovery of motile sperm from thawed human semen. *Archives of Andrology*, 1991; 26: 155-61.

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