

Thermo Cell Transporter 3018

Quality – Made in Germany



Transport incubator for safe and versatile use

Features:

- Factory-aligned temperature at 37°C
- Different temperature possible (to be advised before order)
- Continuous operation with power supply, car voltage or integrated accumulator during transport
- Operation with batteries up to 5 hours*
- Compact design

* at 22°C ambient temperature

Technical data:

- Dimensions: 111 x 111 x 215 (w x d x h in mm)
- Weight: 2.4 kg
- Power supply: 100 – 240 V AC, 50/60 Hz
- Temperature: between 32°C and 42°C available (to be advised before order)
- Accuracy of regulation: $\pm 0.3^\circ\text{C}$
- Operation time:
 - with power supply or 12V car voltage continuously
 - with accumulator approx. 5hrs*
- Heating up time 22°C to 37°C: approx. 1h*
- Charging time of accumulator: approx. 3hrs
- For use with 8 Falcon tubes (round bottom, 13 ml)

Scope of supply:

- Basic unit
- Power supply unit
- Car connection cable

REF 13881



certified company



Rev. 2_11/2013