

Blockthermostat

Quality – Made in Germany



Technical data:

- Dimensions thermostat:
290 mm x 175 mm x 110 mm
(W x D x H)
- Power supply:
100 – 240 V AC, 50/60 Hz, 60 W
- Temperature control between
27°C and 45°C
- Accuracy: +/- 0.2°C
- Temperature distribution: +/- 0.2°C
- Heating up time from 22°C to 37°C:
approx. 40 mins.
- Safety class II

Designed to maintain the temperature of biological materials

Features:

- Temperature adjustable between 27°C and 45°C
- Homogeneous temperature distribution
- Heating up time only approx. 40 min.
- Use of different heating blocks

REF 10719 (without heating block)

REF 13855 (with heating block: REF 10264)

REF 13896 (with heating block: REF 10277)



certified company



Rev. 3_11/2013