

Conical Tube

(#50015/50050/50115/50150/51015/51115/50040/50250)

1. Applications

Conical tubes are widely used in laboratories for multiple purposes, from storage of samples to high speed centrifugation. SPL Life Sciences provides very high strength and non-toxic conical tubes to meet the various demands of scientists.

2. Features

- ❖ Screw cap with tight sealing
- ❖ External graduations with marking area
- ❖ For high speed centrifugation (Cat. No. 50040)
- ❖ Skirted, self-standing type (Cat. No. 50250)
- ❖ Clear polypropylene (Cat. No. 50015, 50050, 50115, 50150, 50040, 50250)
- ❖ Clear polystyrene (Cat. No. 51015, 51115)
- ❖ Provided with plastic rack (Cat. No. 50115, 51115, 50150)
- ❖ For high speed centrifugation (Cat. No. 50040)
- ❖ Sterilized
- ❖ Non - pyrogenic
- ❖ Non - cytotoxic

3. Functional Test

Sterility : Passed
 Cytotoxicity : Passed
 Non-Pyrogenic : Passed
 RNase/DNase Free : Passed

4. Product Range

Cat.No.	Material (Body/Cap)	External Dimension (mm)	Total Volume (ml)	RCF Rating	Units per Rack / Units per Sleeve	Units per Case
50015	PP / HDPE	17 x 120	15	9000	50 / Sleeve	500
50050	PP / HDPE	30 x 115	50	10000	25 / Sleeve	500
50115	PP / HDPE	17 x 120	15	9000	25 / Rack	500
50150	PP / HDPE	30 x 115	50	10000	25 / Rack	300
51015	PS / HDPE	17 x 120	15	3000	50 / Sleeve	500
51115	PS / HDPE	17 x 120	15	3000	25 / Rack	500
50040	PP / HDPE	30 x 115	50	20000	10 / Sleeve	100
50250	PP / HDPE	30 x 116	50	3000	25 / Sleeve	500

5. ISO 9000 Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2008 Quality Management System (K-QA-Q031478) in Korea.


6. Specifications & Technical Drawing

(1) 15 ml Conical Tube (# 50015/50115/51015/51115)

- Specifications

	Dimensions	Material (Resin)	Packaging
Cap	Diameter: 21.30 mm Height: 11.00 mm	High-density polyethylene (HDPE)	<ul style="list-style-type: none"> ▪ Pieces/sleeve: 25 or 50 ▪ Pieces/case: 500
Body	Diameter: 16.90 mm Height: 119.20 mm Volume: 15.00 ml	Polypropylene(PP) or Polystyrene(PS)	

- Technical Drawing




(2) 50 ml Conical Tube (# 50050/50150/50040/50250)

- Specifications

	Dimensions	Material (Resin)	Packaging
Cap	Diameter: 34.00 mm Height: 13.00 mm	High-density polyethylene (HDPE)	<ul style="list-style-type: none"> ▪ Pieces/sleeve: 15 or 25 ▪ Pieces/case: 100 or 500
Body	Diameter: 29.30 mm Height: 113.90 mm Volume: 50.00 ml	Polypropylene(PP) or Polystyrene(PS)	

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr