

PCR Tube Rack

(#80096)

1. Applications

PCR Tube Racks are useful for handling and storage of 0.2 ml PCR tubes

2. Features

- ❖ Autoclavable
- ❖ 5 colors: White, Blue, Green, Red, Purple.


3. Product Range

Cat. No.	Material	Description	Units per Case
80096	PP	8 x 12 (96 Holes)	10

4. ISO 9000 Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2008 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs




6. Specifications & Technical Drawing

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 86.00 mm Width: 123.20 mm Height: 26.00 mm	Polypropylene (PP)	▪ Pieces/Case: 10

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr