

Storage Box

(#80100/81081/81100)

1. Applications

Storage Boxes are useful for storing vials in low temperature environments without distortion of storage boxes.

2. Features

- ❖ For 1.5 ml microtubes
- ❖ Temp. range: -70 ~ +140°C
- ❖ 5 colors: Blue, Green, Orange, Pink, Purple (Cat. No. 81100)
- ❖ Hinged lid (except Cat. No. 81081)

3. Product Range

Cat. No.	Material	Description	Feature	Units per Case
80100	PP	10 x 10 (100 Holes)	Natural Color	30
81081	PP	9 x 9 (81 Holes)	Natural Color	40
81100	PP	10 x 10 (100 Holes)	5 colors	30

4. ISO 9000 Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2008 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs




6. Specifications & Technical Drawing

- Specifications


Dimensions		Material (Resin)	Packaging	
Cat. No. 81081	Cat. No. 80100/81100		Cat. No. 80025	Cat. No. 80081
Length: 133.50 mm	Length: 140.50 mm	Polypropylene (PP)	• Pieces/Case: 40	• Pieces/Case: 30
Width: 133.50 mm	Width: 140.50 mm			
Height: 35.20 mm	Height: 30.00 mm			

- Technical Drawing

Cat. No. 81081



Cat. No. 80100/81100



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr